



GE HealthCare

# Ultrasound

Imaging Accessories & Supplies



Reliable



Compatible



One-Stop Shop

# How to Order



Available 24/7, GE HealthCare Service Shop makes ordering supplies online fast and easy.

- Search real-time inventory of compatible, OEM products
- Take advantage of Auto Reorder for parts you buy regularly
- Batch order multiple items at once with Quick Order
- Get instant order and tracking confirmations
- Pay by PO or credit card or integrate your purchasing system

Visit GE HealthCare Service Shop to register.

1. Visit GE HealthCare Service Shop, navigating to the registration link in the upper left-hand corner of the page.
2. Give us your account details. Pro tip: Provide your customer account number for expedited registration.
3. A welcome message will come to you through the email you provide. Start shopping!

Register now at [ServiceShop.GEHealthcare.com](http://ServiceShop.GEHealthcare.com)

## Other ways to order:



Phone: (866) 281-7545 Options 1, 2, 1, 2

Fax: (877) 279-6990

Email: [gehcaccessorysales@ge.com](mailto:gehcaccessorysales@ge.com)

EDI: Free – contact us for details

Hours: 7AM – 6PM (CST)



# Table of Contents

- Transducer Storage ..... U/S 4 ▶
- ABUS System Consumables ..... U/S 5 ▶
- General Purpose Transducer Covers..... U/S 6 ▶
- Endocavity Transducer Covers..... U/S 7 ▶
- Ultrasound Lotion, Disinfecting Wipes & Gels ..... U/S 9 ▶
- Endocavity Needle Guides..... U/S 12 ▶
- Furniture..... U/S 14 ▶
- Printer Paper & Accessories..... U/S 15 ▶

# Disposable Probe Storage Covers

## Protek™ Non-sterile Probe Covers

- Fits endocavity and general purpose probes
- Helps protect entire probe, including handle
- Designated space for applying disinfection label
- Latex-free



E8013AZ



E8013AY

Cat#	Description	Packaging	Warranty
E8013AZ	Protek™ Non-sterile Disposable Storage Cover – 4" x 16" clear storage cover	100/box 4 boxes/case	H
E8013AY	Protek™ Non-sterile PullUp™ Disposable Storage Cover – 4" x 16" PullUp Storage cover with Secure Tape	200/box	H

## • Ultrasound Probe Rack

Our Ultrasound Probe Rack has been designed to hold a variety of probes.

- Cable support on back plate can hold small, narrow probes ie: MSK
- Includes 4 pre-drilled holes for wall mounting
- Weight – 13.5 lbs
- Measurements – 20 in (L) x 15 in (H) x 5.5 in (D)



E8363JF



Cat#	Description	Warranty
E8363JF	Ultrasound Probe Rack	H

**Compatibility:** Long intracavity, large 4-D OB and standard size probes

# ABUS System Consumables

## ABUS Membranes

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables better image quality during the scan. The membrane is single-use-only which eliminates the risk of cross-contamination and the need for sterilization.

### ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS and somo•v™ ABUS Platinum

Catalog #E8340AB

Curved membrane dimensions:

229 mm (L) x 213 mm (W) x 19 mm (D)

### Polysonic® Ultrasound Lotion and Protex™ Ultrasound Disinfectant Wipes

See details on page U/S 18



E8340AB



E8365EG

E8365EB



E8340AB



E8365EF

E8365EC



E8365EH

GE HealthCare offers an array of consumable supplies for the ABUS (Automated Breast Ultrasound) portfolio of systems including scan membranes, Polysonic® ultrasound lotions, and Protex™ disinfectants.

ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo•v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile)

Compatible for use with either ABUS membrane

Compatible for use with ABUS system

Cat#	Description	Quantity per Box	Warranty
E8340AB	ABUS Single-Use Stabilization Membrane for Invenia™ ABUS 2.0, Invenia™ ABUS, and somo•v™ ABUS (Automated Breast Ultrasound) Platinum systems (non-Sterile)	50/box	W
E8365EC	Polysonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles	12 bottles/box	H
E8365EF	Polysonic Ultrasound Lotion (Non-Sterile) – 1 US gallon with refillable dispenser	4 gallons/case	H
E8365EB	Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5 in x 6.0 in	12 pack/box	H
E8365EG	Protex Ultra Disinfectant Wipes – 75 ct canister pack – 7 in x 9.5 in	8 pack/case	H
E8365EH	Protex Disinfectant Wipes - 80 ct Softpack – 7 in x 10 in	12/case	H

# General Purpose Transducer Covers

## PullUp™ Sterile Latex-Free General Purpose Transducer Covers

- Minimizes cross-contamination
- Protects equipment and reduces post-procedure clean-up
- PullUp™ allows cover to be easily propped open for convenient placement of gel and probe
- PullUp™ handle allows for quick application while maintaining a sterile field
- Made of soft polyurethane
- 3D end for seamless scanning
- Sterile, intended for single-use
- Available in lengths of 24 and 48 in
- Packaged 24 Covers per box, sold per box



E8013AD/AV

Cat#	Description	Packaging	Warranty
E8013AD	5.5 in x 48 in (14 cm x 122 cm) PullUp™ with 3D end, sterile gel and bands	24/box	K
E8013AV	5.5 in x 24 in (14 cm x 61 cm) PullUp™ with 3D end, sterile gel and bands	24/box	K

# Intraoperative Transducer Covers

## Protek™ Medical Sterile Latex-Free Polyethylene Intraoperative Transducer Covers

- Minimizes cross-contamination
- Protects and extends equipment life
- Reduces post-procedure clean-up
- Accordion-folded for easy application
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use
- Intraoperative transducer cover for a sterile environment
- Extends up the length of the transducer cord, allowing up to 8 ft (2.4 m) of sterile protection
- Catalog #E8013AL includes sterile gel packets and bands
- Packaged 12 Covers per box, sold per box



E8013AF



E8013AL

Cat#	Description	Packaging	Warranty
E8013AF	6 in tapered to 3 in x 96 in (15 cm-8 cm x 244 cm) Probe cover and bands	12/box	K
E8013AL	6 in tapered to 3 in x 96 in (15 cm-8 cm x 244 cm) Probe cover, boot, sterile gel, and bands	12/box	K










# Endocavity Transducer Covers – Non-Sterile


## Endocavity Probe Covers – Non-Sterile

- Used for OB/GYN scanning procedures
- Helps to protect against cross-contamination
- Available for Auto-Reorder



	Cat#	Latex Free	Description	Bands Included	Gel Included	Packaging
	E8013AK	X	12 in (30.48 cm)	X		50/bag
	E8013AR		0.8 x 7.9 in (2 x 20 cm)	X		24/box
	E8013AM		1 x 7.9 in (2.6 x 20 cm)	X		24/box
	E8013AE		1.4 x 7.9 in (3.5 x 20 cm)	X		24/box
	E8013AN		1.4 x 7.9 in (3.5 x 20 cm)			100/bag
	E8013AP		1.4 x 7.9 in (3.5 x 20 cm)			100/bag; 5 bags/box
	E8386CG	X	2.37 in tapering to 1.75 in (W) x 9.5 in (L) (60 mm tapering to 44 mm (W) x 241 mm (L))		X	100/bag

### Endocavity Probe Cover Bands


	Cat#	Latex Free	Description	Packaging
	E8013AG	X	.75 in (1.9 CM)	500/bag


NOTE: All Endocavity Probe Covers have Warranty Code K.

# Endocavity Transducer Covers – Sterile

## Protek™ Endocavity Probe Covers – Sterile

- Used for OB/GYN scanning procedures
- Helps to protect against cross-contamination
- Available for Auto-Reorder



	Cat#	Sterile	Latex Free	Description	Bands Included	Gel Included	Packaging
	E8013AH	X	X	5 in tapered to 2 in x 24 in (13cm-5cm x 61cm)	X	X	24/box



# Ultrasound Lotion

## Parker Labs Polysonic® Ultrasound Lotion

POLYSONIC Ultrasound Lotion is an acoustic coupling lotion that is hypo-allergenic and has a light pleasing scent.



E8365EF

E8365EC

Cat#	Description	Packaging	Warranty
E8365EC	PolySonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles	12/box	H
E8365EF	PolySonic Ultrasound Lotion (Non-Sterile) – 1 US gallon	4 gal/case	H

**Compatibility:** Used with GE ABUS system

# Ultrasound Disinfectant Wipes

## Parker Labs Protex Ultrasound Disinfectant Wipes

Protex™ one-step cleaner disinfectant kills more than 30 pathogens including HIV, H1N1 and MRSA, while remaining safe for even the most sensitive surfaces. The pre-moistened wipes are easy and convenient to use.



E8365EG

E8365EB



E8365EH

Cat#	Description	Packaging	Warranty
E8365EB	Protex Ultra Disinfectant Wipes – 60 ct Soft-pack w/stay-put adhesive backing – 6.5 in x 6.0 in	12 pack/box	H
E8365EG	Protex Ultra Disinfectant Wipes – 75 ct canister pack – 7 in x 9.5 in	8 pack/case	H
E8365EH	Protex Disinfectant Wipes - 80 ct Softpack – 7 in x 10 in	12/case	H

**Compatibility:** Visit [gehealthcare.com/ultrasound](http://gehealthcare.com/ultrasound) for detailed information.

# Ultrasound Gels

## Parker Labs Aquasonic 100 Gel

- Original gel for diagnostic and therapeutic applications
- Tested and approved for use with GE transducers
- Blue colored aqueous gel formula is non-staining, non-sensitizing, and non-irritating
- Hypoallergenic, bacteriostatic
- No salt or formaldehyde; not a spermicide
- 5 liter jug with dispenser, 250 ml bottles, or 20 g sterile overwrapped foil pouches



E8365BA

E8365AF

E8365BJ

Cat#	Description	Packaging	Warranty
E8365BA	Aquasonic 100 Gel (250 ml bottles)	12/case	H
E8365AF	Aquasonic 100 Gel (5 L jug – includes dispenser)	1/case	H
E8365BJ	Aquasonic 100 Gel (20 g pouches – single-use)	48/box	H

## Parker Labs Aquasonic Clear Gel

- Same formula as original Aquasonic 100 Gel, but without the perfume or color
- Clear, non-scented aqueous, water-soluble gel formula is non-staining, non-sensitizing and non-irritating
- Hypoallergenic, bacteriostatic
- Offered in 250 ml bottles or 5 liter jug with dispenser



E8365CB

E8365CA

Cat#	Description	Packaging	Warranty
E8365CA	Aquasonic Clear Gel (250 ml bottles)	12/case	H
E8365CB	Aquasonic Clear Gel (5 L jug – includes dispenser)	1/case	H

**Parker Labs SCAN<sup>®</sup> Gel**

- Dual viscosity formula allows SCAN gel to be applied as a gel and scan as a fluid
- Less viscous than Aquasonic 100 for ultrasound procedures where a less viscous gel is desired
- Green colored gel is completely water-soluble, will not stain clothing, damage equipment, or irritate skin
- Hypoallergenic, bacteriostatic
- Offered in 250 ml bottles, 1 gallon jug with dispenser



E8365BC

E8365BD/BK

Cat#	Description	Packaging	Warranty
E8365BC	SCAN Ultrasound Gel (8 oz. bottles)	12/case	H
E8365BD	SCAN Ultrasound Gel (1 gal. jug)	1/case	H

## Standoff Pads

**Parker Labs Aquaflex<sup>®</sup> Standoff Pads**

- Aqueous, flexible, disposable ultrasound standoff pad for use in difficult to visualize and near field areas
- Requires no gel, except in prolonged procedures
- Cushions sensitive areas for comfortable scanning
- Hypoallergenic, bacteriostatic
- Pads are disposable, which helps to prevent cross-contamination
- Each pad measures 2 × 9 cm
- Packaged 6 pads per box, sold per box



E8317C

Cat#	Description	Packaging	Warranty
E8317C	Aquaflex Standoff Pads – 2 × 9 cm	6/box	H

# Endocavity Needle Guides

## Protek™ Disposable Endocavity Needle Guide (Single Piece)

- For precise needle placement during ultrasound-guided procedures
- Directs instruments according to on-screen biopsy guidelines
- Easily snaps onto probe
- Disposable design of guide helps to reduce the risk associated with cross-contamination
- Smooth profile ensures patient comfort
- Wide funnel for optimal needle entry
- Requires no disinfecting
- Fits needle gauges 16 or thinner
- Sterile, intended for single-use
- 24 disposable procedure kits



E8013AS



E8013AT

Cat#	Description	Packaging	Warranty
E8013AS	E8C/E721/E8C-RS/IC5-9H/MTZ Sterile 16Ga single-use needle guide Includes 24 disposable biopsy needle guides	24/box	K
E8013AT	E8C/E721/E8C-RS/IC5-9H/MTZ Biopsy Kit Sterile 16Ga single-use needle guide, 3.5cm x 20cm latex cover, sterile gel, and bands	24/box	K

**Compatibility:** GE HealthCare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ Transducers

**Protek™ Disposable Endocavity Needle Guide**

- For precise needle placement during ultrasound-guided procedures
- Easily snaps onto probe
- Disposable design of guide reduces risks associated with cross-contamination
- Smooth profile ensures patient comfort
- Wide funnel for optimal needle entry
- Requires no disinfecting
- Fits needle gauges 16 or thinner
- Sterile, intended for single-use
- 25 disposable procedure kits




E8013AQ

Cat#	Description	Packaging	Warranty
E8013AQ	Sterile 16Ga single-use needle guide	25/box	K

**Compatibility:** GE HealthCare 618E, TR7 Transducers

**Protek™ Disposable Endocavity Needle Guides**

- For precise needle placement during ultrasound-guided procedures
- Directs instruments according to on-screen biopsy guidelines
- Easily snaps onto probe
- Smooth profile ensures patient comfort
- Wide funnel for optimal needle entry
- Requires no disinfecting
- Fits needle gauges 16 or thinner
- Sterile, intended for single-use
- Each Disposable Biopsy Needle Guide Starter Kit (E8013AW) includes 24 sterile kits:
  - Sterile 16Ga needle guide, 5 in tapered to 2 in x 24 in soft cover, 2.0cm x 20cm latex cover, and bands



E8013AW

**Protek™ Reusable Endocavity Needle Guide**

- For precise needle placement during ultrasound-guided procedures
- Directs instruments according to on-screen biopsy guideline
- Attaches securely to probe with ease (style varies)
- High-quality, stainless steel
- Safe for sterilization (steam sterilize/autoclave)
- Non-sterile, reusable



E8013AX

Cat#	Description	Packaging	Warranty
E8013AW	Sterile 16Ga needle guide, 5 in tapered to 2 in x 24 in soft cover, 2.0cm x 20cm latex cover, and bandst	24/box	H
E8013AX	Non-sterile 16Ga stainless steel needle guide (reusable)	1	H

**Compatibility:** GE HealthCare BE9C, BE9CSTransducers

# Ultrasound Furniture

## GE Anti-Fatigue Mat


- Helps to reduce body stress and fatigue resulting from standing for long periods
- One of the most comfortable, buoyant, anti-fatigue floor mats on the market
- Durable construction; sealed to prevent moisture absorption and facilitate cleanup



E7058A



E7058AB



E7058AC

Cat#	Description	Warranty
E7058A	GE Anti-Fatigue Floor Mat (Blue 3' x 5' x 5/8") – 20 lb	H
E7058AB	GE Anti-Fatigue Floor Mat (Gray 3' x 5' x 5/8") – 20 lb	H
E7058AC	GE Anti-Fatigue Floor Mat (Blue 2' x 3' x 5/8") – 9 lb	H

## Physician Mobile Armchair

- Padded arms
- Charcoal gray
- Adjusts from 17.7 to 21.7 in (45 to 55 cm)
- Weight 31 lb (14 kg)
- Chair contains some ferro-magnetic material
- Not safe for use in an MRI Suite's Zone IV, use in control room only



E80331BE

Cat#	Description	Warranty
E80331BE	Physician Mobile Armchair	H

# Sony Printer Media Supplies

## Sony Printer Media and Supplies

• GE HealthCare offers a variety of high-quality Sony printer papers and supplies to meet all of your printing needs. Please refer to the table below to find the Sony media recommended for your printer



E8310JT



E8310JS



E8310KP



E8310KS



E8310LA



E8310LB



E8310LC



E8310GA

Cat#	Cross-Ref	Description	Compatible with Sony Printers	Prints Per Pack or Roll	Sold Per	Warranty
E8310JT	UPP-210HD	Sony High-Density Black & White Media (130 prints per roll; 1 case contains 5 rolls)	UP-991AD, UP-971AD, UP-910, UP-930, UP-960, UP-980	130	1 case = 5 rolls	H
E8310JS	UPP-210SE	Sony Standard Black & White Media (130 prints per roll; 1 case contains 5 rolls)	UP-991AD, UP-971AD, UP-910, UP-930, UP-960, UP-980	130	1 case = 5 rolls	H
E8310KP	UPC-21L	Sony A6 Color Print Pack (200 prints of paper and ribbon – 50 sheets x 4 packs)	UP-D25MD, UP-25MD, UP-21MD, UP-D23MD, UP-20	200	1 carton = 6 packs	H
E8310KS	UPC-55	Sony A5 Color Print Pack (2 individual packs of 100 sheets each and 2 ribbons)	UP-D55, UP-D55MD	200	1 carton = 5 packs	H
E8310LA	UPP-110HA	Sony High-Density Black & White Media (240 prints per roll; 1 case contains 10 rolls)	UP-880, UP-890MD, UP-D890	240	1 case = 10 rolls	H
E8310LB	UPP-110HD	Sony High-Density Black & White Media (270 prints per roll; 1 case contains 10 rolls)	UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	270	1 case = 10 rolls	H
E8310LC	UPP-110HG	Sony High-Glossy Black & White Media (240 prints per roll; 1 case contains 10 rolls)	UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	240	1 case = 10 rolls	H
E8310GA	UPP-110S	Sony Standard Black & White Media (270 prints per roll; 1 case contains 10 rolls)	UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897MD, UP-D898MD, UP-X898MD	270	1 case = 10 rolls	H



## Specifications

- All specifications are subject to change without notice
- The product and brand names herein are the trademarks or registered trademarks of their respective owners

## Pricing

- Please call to obtain current pricing. All pricing and product availability is subject to change without notice
- Price invoiced will be the price in effect at the time of order entry
- Quotations available upon request

## Blanket Purchase Orders

You can save time and money by ordering your consumable products with our Blanket Order program. Blanket orders must be scheduled for a year or more and a signed purchase order agreement must be on file. Contact us for additional details.

## Restocking Fee

- Returns are subject to, but not limited to, a minimum 15% restocking fee
- Certain items, such as sterile and environmentally controlled products within this catalog, are non-returnable and nonrefundable unless defective and have been noted as such

## Product Returns

- Immediately upon receipt, please inspect all deliveries for order accuracy and condition
- For products shipped in error, please obtain return authorization from the Accessories and Supplies Call Center within thirty (30) days of delivery
- Power Quality items are nonreturnable and nonrefundable

## Contact Us for General Information and Warranty Terms

### In U.S.

GE HealthCare  
Accessories and Supplies Customer Support Center  
P.O. Box 414, W-544  
Milwaukee, WI 53202

Phone: (866) 281-7545 – opt 1, opt 2, opt 1, opt 2

Fax: (877) 279-6990

Online: [www.gehealthcare.com/accessories](http://www.gehealthcare.com/accessories)


Email: [gehaccessorysales@ge.com](mailto:gehaccessorysales@ge.com)

### All other regions

Contact your local GE Representative.

# Save time. Order accessories online with Service Shop.

## Auto-Reorder

 Auto-Reorder

**Set quantity, schedule frequency**

Every 3 Months ▾

Customize the quantity and order frequency that works for you.

 Add to Cart

**Add to cart**

It's that easy. Your first order has been placed. The product will be auto-delivered at your preferred delivery frequency.



**Manage reorders**

In advance of each delivery, we will send you a reminder email for your upcoming delivery. Skip your deliveries or cancel your subscriptions at any time.



## Quick Order

Know what you need to order? Batch order multiple items at once with Quick Order on Service Shop. Enter a list of item numbers or paste from an excel spreadsheet, select quantity and add to cart. Spend less time ordering!



## Payment Options

Pay by PO or credit card or integrate your purchasing system.

## Register now at [ServiceShop.GEHealthcare.com](https://ServiceShop.GEHealthcare.com)

1. Visit GE HealthCare Service Shop, navigating to the registration link in the upper left-hand corner of the page.
2. Give us your account details. Pro tip: Provide your customer account number for expedited registration.
3. A welcome message will come to you through the email you provide. Start shopping!

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator dedicated to providing integrated solutions to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier.

Learn more at [www.gehealthcare.com](http://www.gehealthcare.com).

GE HealthCare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE HealthCare representative for the most current information and availability. Please visit [www.gehealthcare.com](http://www.gehealthcare.com). Data subject to change.

©2023 General Electric Company.

GE and the GE Monogram are trademarks of General Electric Company.

GE is a trademark of General Electric Company used under trademark license.

All other product names and logos are trademarks or registered trademarks of their respective companies.

JB00640US Rev. 09/23



GE HealthCare