



Sustainable magnetic resonance imaging solutions for a resilient tomorrow

SIGNA™ Premier/SIGNA™ Premier Evo





Creating a more sustainable future requires us to care for the planet and its inhabitants.

It is essential that we continue to drive progress toward early, precise, and accessible diagnosis and treatment of more patients. For the planet, it is critical that we do so with a reduced impact on precious and rare resources that are imperative to life. We believe that the advancement of precision health, greater digitization of healthcare, and increased access to quality care are fundamental to accomplishing this goal.

We support carbon policies that reduce greenhouse gas emissions and promote sustainable development. We are committed to achieving net zero by 2050 and are part of the UN-backed “Race to Zero,” with a goal of reducing emissions based on the Paris Agreement. We’ve also set a public goal to achieve a 50% reduction in our own operational emissions by 2030. As a result of these efforts, we want to enable a more sustainable health system by addressing not only the environmental impacts of our products but also the challenges healthcare professionals and their patients face with resilient, digital options.



We are committed to achieving **net zero** emissions by 2050.

We’ve set a public goal of a **50% reduction** in our own operational emissions by 2030.



We deliver sustainable, intelligently efficient solutions for a resilient tomorrow.

Building a healthier world to help improve access to care and enable better patient outcomes.



Green

Using fewer resources for a healthier planet.



Digital

Transforming healthcare through innovation.



Resilience

Building flexibility and dependability across healthcare systems.



SIGNA™ Premier/SIGNA™ Premier Evo helps create a resilient tomorrow.

Our 3.0T MR system, SIGNA™ Premier/SIGNA™ Premier Evo, and its services help ensure that radiology professionals and the patients they serve have the technology necessary to create a sustainable and resilient tomorrow.

Designed for research and beyond, SIGNA™ Premier delivers an imaging experience that helps differentiate services through technology, clinical capability, and patient comfort. SIGNA™ Premier Evo is a simple upgrade that widens the clinical capabilities of legacy MR systems, extends the legacy system's life expectancy, and reduces its carbon footprint. Keeping the current magnet eliminates waste, reduces helium consumption, and lowers reconstruction costs.

Reducing environmental impact

- The life expectancy of a GE HealthCare 3.0T magnet is up to 30 years.
- GE HealthCare MR systems are eligible for refurbishment, reuse, or recycling at the end of their product life.
- 30% less helium is required, and there is no helium refill onsite during installation.¹

Improving outcomes

- Human Connectome research performance with SuperG Gradient enables exceptional gradient stability for the most demanding fMRI exams.
- The AIR™ family of products delivers clinical versatility and comfort, intelligent productivity improvements, and consistently extraordinary image quality.
- Operational excellence is maintained using digital tools and subscription models to ensure the system stays up to date.



¹ Compared to previous generation MR system.



Contributing to a healthier planet

More than half of the healthcare sector’s climate footprint, approximately 53%, is attributable to energy use.² As a result, we have strengthened our commitment to environmentally conscious design and sustainable practices across our product manufacturing, sourcing, distribution, installation, and service operations. This includes improving energy efficiency, optimizing the use of limited or rare materials, providing digitally enabled and remote predictive and maintenance service throughout the product lifespan, and offering refurbishment and recycling options at the end of product life.

GE HealthCare environmental management system is ISO 14001 certified.

Our production and service operations align to ISO 14001 standards.

We’re committed to environmental product design.

This product conforms with IEC60601-1-9:2007.

Materials

GE HealthCare reviews the environmental aspects of the material supply used within our products to increase recyclability and decrease the use of hazardous substances, when possible.

Recyclable

84% of the MR magnets returned are recycled into raw materials.³

Our MR system covers are completely recyclable plastics.

Reduce the use of hazardous substances

REACH (EC) 1907–2006

Manufacturing

Through our environmental reviews, we also focus on implementing renewable energy and reducing waste.

Reducing electricity

The system uses Ultra High Efficiency (UHE) technology to generate exceptional imaging performance using about half the amount of electrical current than a contemporary system.

² Health care climate footprint report | Health Care Without Harm (noharm-uscanada.org)

³ Data on file.



Product utilization

Our imaging products are designed to help enable energy efficiency through dedicated features and advanced applications to reduce the environmental impact.

Ergonomically designed

Patient setup and positioning

Transforming from a 60 cm to a 70 cm bore widens the cross-section by 36%, giving patients extra room during their scan.

Reduce claustrophobia rejection rates by 90% while improving general comfort and patient satisfaction.⁴

AIR™ Recon DL enables shorter scan times, reducing the time spent on the table.

AIR™ Coils are designed to be flexible and 50% lighter to provide better coverage than traditional coils. They allow greater patient positioning freedom and patient comfort, as well as reduced burden on the technologists lifting the coils.⁵


AIR Touch™ offers 59% less patient setup time and 37% less table time for the patient.⁶

AIR x™ automated MR slice prescription reduces setup time and provides reproducible planning to ensure exam consistency.

Reduce staff burden

Increase productivity and streamline workflows with shorter scan times.

Optional detachable eXpress Patient Table improves patient handling and comfort, allowing for faster patient setup.



AIR™ Anterior Array (AA) Coil

⁴ Claustrophobia rate comparison head-first/feet-first.

⁵ AIR™ Coils data on file (2018).

⁶ Compared to conventional technology. Data on file.



Packaging

GE HealthCare imaging equipment has a robust and multi-sourced supply chain for systems and spare parts across all product portfolios.

Improved packaging

We are replacing our wood and corrugated cardboard packaging with paper, increasing the amount of recyclable packaging.

Product utilization (Cont.)

Reduce noise

Silenz is a 3D Zero-TE sequence comprised of high-bandwidth excitation and reduced gradient-switching radial acquisition that drastically reduce noise level from an ear splitting, motorcycle-level 91 dB to within 3 dB of scan room ambient sound. In addition, Silenz has added flexibility in sequence prescription to enable faster scan times.

Guidance for product utilization

Instructions are provided for use of the equipment to minimize the environmental impact during installation, use, and operation.

Reduce energy consumption during use

Utilize standby power mode to reduce energy consumption by 35% when the system is idle.⁷

AIR Touch™ offers 59% less patient setup time and 37% less table time for the patient.⁸

Power consumption⁹

Off mode: 8.2 kW
Standby (no scan): 25.8 kW
Scan mode: 41.8 kW

Reduce consumable energy utilization

1115 L of helium are recovered per system.¹⁰

Our helium recovery systems reduce any escaping gas during cool-down, ramp boil off, and training quench.

⁷ Compared to energy consumption when the system is in scan mode.

⁸ Compared to conventional technology. Data on file.

⁹ Per COCIR Self-regulatory initiative for medical imaging equipment, over a 24-hour period, with 12 hours of active day and 12 hours of inactive night scenario.

¹⁰ Value based upon factory recovery system. Data on file.



End of product life

We are increasingly putting our retired products' materials back into the supply chain to maximize efficient use and minimize unnecessary waste. This circularity model enables our imaging products to extend their clinical impact through longer lifespans while reducing the environmental footprint. Additionally, we offer our customers partnered support for upgrades and services throughout a product's lifespan to maintain optimal performance and help drive better patient outcomes.

Our refurbishment programs involve an extensive inspection and testing process, designed to bring equipment back to its original certified manufacturing specifications. If the system is not suitable for refurbishment, eligible parts are harvested for reuse after quality and performance testing, while the rest are returned to dedicated recycling facilities.

Product utilization (Cont.)

Guidance for end of lifecycle

Equipment instructions are provided to minimize the environmental impact for disposal or recycling.

Upgradeable hardware and software options are provided as a solution to extend the product lifespan.

SIGNA™ Premier Evo widens the bore on legacy MR systems to 70 cm, accommodating more patients and extending the life of the system.

Gain a 40% faster ROI vs. a replacement with a SIGNA™ Premier Evo transformation, and reduce your carbon footprint by keeping your current advanced GE HealthCare MR magnet.¹¹

Parts harvesting and refurbishment options are provided to reduce waste and environmental impacts while extending imaging access to less advantaged regions.

MR system parts are eligible for assessment through the refurbishment program, in which they are assessed for refurbishment, harvesting, or recycling at the appropriate time in the lifespan.¹²

94–96% of most systems are reused, refurbished, or recycled, extending the lifetime of each product.¹²

Waste reduction

This system is in accordance with Waste Electrical and Electronic Equipment (WEEE) regulations.

¹¹ Based on historical ROI of full system upgrades vs full replacement.

¹² 100% of products within MR, CT, nuclear medicine, PET, and X-ray are eligible for refurbishment, although whether a system is actually refurbished versus harvested for parts or otherwise recycled or reused is dependent on the state of the system when GE HealthCare takes possession of it. Data on file.



Digitizing healthcare through transformative innovations for a resilient tomorrow

We are committed to investing in digital capabilities that help accelerate clinical decision making, optimize imaging operations, and drive efficiencies in exam workflows, all of which can improve patient outcomes. Enabling digital transformation will further enhance our predictive and maintenance service operations for the life of your products.

We are also dedicated to driving a more resilient and sustainable future in healthcare. Many factors, including the pandemic, climate-related weather disasters, and supply-chain issues amplified this need. Managing operations through these challenges requires resilience and perseverance.

Advancing clinical outcomes

Advanced applications and cutting-edge AI tools provide personalized data to drive actionable insights, helping healthcare professionals make fast, accurate clinical decisions for care pathways.

Gain actionable clinical insights quicker for earlier diagnosis

Scan up to 50% faster and reduce backlogs with AIR™ Recon DL.¹³

Increase SNR and reduce scan time in a single exam.

Offer greater consistency and stability for advanced fMRI and other demanding applications.

Keep your imaging equipment up to date with advanced clinical applications

Smart Subscription protects equipment from obsolescence and keeps the system at its best. It improves patient outcomes and productivity due to improved functionality and easy access to innovation.

Upgrade your clinical applications with SIGNA™ Continuum for all existing installed SIGNA™ MR systems.

¹³ Compared to conventional technology. Data on file.



Advancing clinical outcomes (Cont.)

Help improve patient outcomes with improved image quality

80% of cases see improved SNR without added scan time with AIR™ Recon DL.¹⁴

Using a deep learning-based reconstruction algorithm, AIR™ Recon DL improves SNR by making use of the raw data to remove image noise and ringing.

Drive advancements with precision health

Wing-to-wing clinical solutions from setup to report help reduce scan times and increase precision health, including anatomy-dedicated post-processing tools and quantitative tools for measuring and assisting diagnosis.

¹⁴ Results may vary.



Optimizing imaging operations

Our AI-based and advanced digital solutions are designed to increase efficiencies across the radiology spectrum without increasing the administrative and training burden on radiologists and technologists.

Increase productivity and consistency

Reduce scan time up to 50% with AIR™ Recon DL.¹⁵

AIR x™ automated MR slice prescription reduces setup time and provides reproducible planning to ensure exam consistency.

GE HealthCare's Imaging Protocol Manager remotely updates and manages protocols between facilities on MR systems to help deliver consistent image quality and optimal patient care.

Gain data intelligence and actionable insights across the radiology department to increase productivity with Imaging Insights.

MR Performance Excellence can help increase capacity by optimizing equipment usage.

Leverage on-demand or scheduled virtual clinical applications training with GE specialists to support staff, enabled by Digital Expert Access.

¹⁵ Compared to conventional technology. Data on file.



Optimizing imaging operations (Cont.)

Reduce downtime

The OnWatch™ remote monitoring system reduces unplanned down time by 35% and service calls by 39%, with an 11% reduction in on-site repair time.¹⁶

This, in turn, helps reduce travel and carbon footprint, as well as overall energy and waste, by keeping our systems optimally running.

Improved thermal performance allows the magnet to stay cold longer during extended outages, and the system can immediately return to scanning when the power goes back on.

Backup UPS for power outages enable exam completion, while long-term UPS protects the system from helium loss. Backup generator cryogenic management allows the system to resume operations when feasible.

Cybersecurity

SIGNA™ Premier and SIGNA™ Premier Evo incorporate a broad assortment of capabilities to enable privacy and security. The layered approach of defense in depth limits the risk that the failure of a single security safeguard will allow compromise of the system. GE HealthCare's Design Engineering Privacy and Security (DEPS) process follows GDPR, HIPAA, NIST 800-53, NIST 800-30, ISO 27001, and NIST CSF requirements.

¹⁶ Compared to previous generation software. Data on file.



Enabling intelligent exam workflows

Intelligent automation features help drive consistency, enable fast, easy exams, and improve workflow with fewer resources, all while achieving similar or improved outcomes.

Reduce setup time

AIR Touch™ reduces setup time by 59%.¹⁷

AIR x™ allows five times faster set up with four times fewer mouse clicks.¹⁸

Reduce exam time

AIR Touch™ reduces patient table time by 37%.¹⁷

Reduce scan time up to 50% reduction per patient with AIR™ Recon DL, enabling improved workflow and efficiency.¹⁹

Reduce overall scan times without compromising image quality with HyperSense, which can be used in 88% of all clinical procedures.

Ease of use

AIR™ Coils optimize SNR, enabling positioning for any patient shape, while anatomy-dedicated post-processing tools unique to the specific organ system bring precision medicine into the MR suite.

Accelerate emergency egress with the eXpress Patient Table, which can be undocked and removed by one user in under 30 seconds, typically.

Cleanability

Our equipment is designed to be cleaned and disinfected easily. We continue to test and approve new cleaning and disinfecting agents. Visit [Cleaning.GEHealthCare.com](https://www.gehealthcare.com/cleaning) for updates.

¹⁷ Compared to previous generation software. Data on file.

¹⁸ Comparison of automated workflow with AIR x™ versus traditional setup process. Data on file.

¹⁹ Compared to conventional technology. Data on file.



Building a healthy world to help enable better patient outcomes.

GE HealthCare is a member of COCIR, the European Trade Association representing the medical imaging, radiotherapy, health ICT, and electromedical industries.²⁰

²⁰<https://www.cocir.org/about-cocir/members.html>

Not all products or features are available in all geographies. Check with your local GE HealthCare representative for availability in your country. Not all features are included in the standard system configuration. Check with your local GE HealthCare representative.

© 2023 General Electric Company. GE, the GE Monogram, SIGNA, AIR, AIR Touch, OnWatch, and AIR x are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

JB27537XX