

GoldSeal™ Ultrasound

General Imaging family of products



LOGIQ E10s

LOGIQ E10

LOGIQ Fortis



LOGIQ S8



LOGIQ P9 / P10



Vscan Air

As the demand for medical imaging continues to grow, pre-owned systems help meet an important need for those seeking cost-effective imaging solutions. With the GE HealthCare GoldSeal refurbishment program for Ultrasound, you can feel confident you are purchasing a quality system at a great value that only the original equipment manufacturer (OEM) can deliver.

With thousands of GoldSeal ultrasound systems in service around the world, you know our reputation. You can trust us to provide quality imaging equipment so that you can provide quality patient care.

Confidence comes standard with every GoldSeal system. It will look like a new machine and it will perform like a new machine. We guarantee it.

GoldSeal pre-owned, refurbished systems offer you a smart use of resources, as well as:

- **Selective process:** Stringent selection standards ensure that only those systems with well-known and acceptable service histories qualify for GoldSeal.
- **Quality:** Refurbishing completed by OEM factory-trained technicians who ensure all OEM specifications are met.
- **One-year Warranty:** Same warranty as on most new systems; service contracts available.
- **Service/support:** Online assistance with questions, and local service.
- **Up-to-date technology:** Refurbishing includes installation of latest possible software release and original OEM parts.
- **Training:** Operation and application training is available, with optional CE Tech training credits.

Helping meet your imaging needs at an affordable price

LOGIQ™ E10/E10s

LOGIQ E10/E10s is designed to facilitate accurate diagnosis, two-dimensional shear wave elastography quantifies tissue elasticity and provides color-coded elastograms. This advanced tool helps reduce the need for invasive procedures and empowers clinicians to be more confident in their patient management decisions.

LOGIQ Fortis

LOGIQ Fortis is a powerful streamlined ultrasound built to help users deliver on the promise of confident care in multiple clinical settings. Small and sleek enough to move to patients across departments and equipped with our most powerful technology. Exceeding expectations to facilitate next generation imaging LOGIQ Fortis features cSound Architecture, which combines versatile XDclear™ probes, cSound Imageformer and new, advanced Speckled Reduction Imaging (SRI) technology resulting in increase of processing power that delivers enhanced data throughput for exceptional image quality, clarity, and confidence.

LOGIQ S8

LOGIQ S8 with XDclear provides superb image quality, simplified workflow and can be suited to your needs; built upon Agile Acoustic Architecture technology with XDclear transducer technology helps you achieve superb images even on difficult to image patients using innovative tools. Tailor the system to your practice with our selection of sophisticated clinical tools, educational opportunities, service, and financial solutions.

LOGIQ P9/P10

LOGIQ P9/P10 is a highly capable ultrasound system that provides excellent image quality and productivity through easy to use tools across a wide range of applications. Imaging engine of LOGIQ P9/P10 delivers crisp images in a mobile, ergonomically designed console. Automation and personalization capabilities make the system fast and simple, and personalized settings and innovative tools on the system help you quickly get to the root of each patient's problem.

Vscan Air™

Vscan Air is a battery-operated general-purpose diagnostic ultrasound imaging system. It enables ultrasound imaging guidance, visualization and measurement of anatomical structures and fluid. Vscan Air consists of a dual-headed probe, which integrates both curved and linear array transducers, and an app that can be installed on Android™ or iOS® mobile devices. Its pocket-sized portability and simplified user interface enables integration into training sessions and examinations in professional healthcare facilities.



GoldSeal systems are quoted subject to availability. To confirm current availability of systems or additional options in your region, contact your GE HealthCare representative.

© 2023 GE HealthCare.

Not all products or features are available in all markets. Contact a GE HealthCare Representative for more information. Please visit <http://www.gehealthcare.com>. Data subject to change. LOGIQ, Vscan Extend, XDclear, and GoldSeal are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license. Android is a trademark of Google, Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries

JB11523XX

For more information, please contact your GE HealthCare sales representative or visit online at gehealthcare.com/goldseal