

GoldSeal™ Ultrasound

Cardiovascular family of products



Vivid E95



Vivid S70N



Vivid T9



Vivid iq



Vscan Air

The Vivid product line is specifically designed to meet the imaging and workflow needs in cardiovascular ultrasound. Vivid technologies and workflow tools provide fast, accurate visualization and quantification, while reducing clicks and repetitive tasks.

As the demand for medical imaging continues to grow, pre-owned systems help meet an important need for those seeking cost effective imaging solutions. With the GE HealthCare GoldSeal refurbishment program for ultrasound, you can feel confident you are purchasing a quality system at a great value that only the original equipment manufacturer (OEM) can deliver.

With thousands of GoldSeal ultrasound systems in service around the world, you know our reputation. You can trust us to provide quality imaging equipment so that you can provide quality patient care.

Confidence comes standard with every GoldSeal system. It will perform comparably to new equipment.

GoldSeal pre-owned, refurbished systems offer you a smart use of resources, as well as:

- **Selective process:** Stringent selection standards ensure that only those systems with well-known and acceptable service histories qualify for GoldSeal
- **Quality:** Refurbishing completed by OEM factory-trained technicians who ensure all OEM specifications are met
- **One-year warranty:** Same warranty as on most new systems; service contracts available
- **Service/support:** Online assistance with questions, and local service
- **Up-to-date technology:** Refurbishing includes installation of latest possible software release and original OEM parts
- **Training:** Operation and application training is available, with optional CE Tech training credits

Helping meet your imaging needs at an affordable price

Vivid™ E95

Get to the heart of the issue. The premium Vivid E95 quickly brings you quantifiable information across a wide variety of patients. GE HealthCare's propriety cSound™ software beamformer and image reconstruction platform, combined with GE HealthCare's XDclear™ technology, captures data from every channel in the probe. This delivers exceptional image quality with impressive penetration and/or high sensitivity while maintaining high spatial resolution. Count on a full suite of simple, intuitive tools to help make your work easy and efficient.

Vivid S70N

Vivid S70N is a very mobile, robust cardiovascular ultrasound system with 4D capability. cSound software-based beamforming platform delivers excellent image quality in cardiovascular and shared service applications. 4D TEE operation helps broaden utilization of the Vivid S70N across an organization. Quantification with Vivid S70N is all about helping providers evaluate problems and pursue the path forward. Count on a full suite of intuitive tools to make your work easy and efficient. In addition, its small size, weight and battery stand-by capability deliver a highly mobile solution.

Vivid T9

Ergonomically-enhanced and ideally designed for cardiovascular practices and shared services, Vivid T9 combines the established, intelligent cardiac imaging capabilities of GE HealthCare Vivid systems with exceptional shared services performance of the LOGIQ™ systems. A truly hybrid cardiovascular ultrasound system that is rugged, reliable, robust, with exceptional ergonomic design, and a rich feature set.

Vivid iq

GE HealthCare's compact directly addresses your changing needs with a combination of portability and power that fits right in — wherever your journey takes you — adult echo, interventional, ICU/CCU, ER, pediatrics, and office-based practice. Count on Vivid iq's combination of power and portability to help you achieve the speed and quality in echo this environment demands.

Vscan Air™

Vscan Air is a battery-operated general-purpose diagnostic ultrasound imaging system. It enables ultrasound imaging guidance, visualization and measurement of anatomical structures and fluid. Vscan Air consists of a dual-headed probe, which integrates both curved and linear array transducers, and an app that can be installed on Android™ or iOS® mobile devices. Its pocket-sized portability and simplified user interface enables integration into training sessions and examinations in professional healthcare facilities.



GoldSeal systems are quoted subject to availability. To confirm current availability of systems or additional options in your region, contact your GE HealthCare representative.

© 2023 GE HealthCare.

Not all products or features are available in all markets. Contact a GE HealthCare Representative for more information. Please visit <http://www.gehealthcare.com>. Data subject to change. cSound, GoldSeal, LOGIQ, Vivid, Vscan Extend, and XDclear are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license. Android is a trademark of Google, Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries.

JB03081XX

For more information, please contact your GE HealthCare sales representative or visit online at gehealthcare.com/goldseal