

Key Versana Features



	Versana Premier™* Powerful. Versatile. Productive.	Versana Balance™* Competent. Efficient. Practical.	Versana Essential™* Simple. Effective. Accessible.	Versana Active™* Advanced. Capable. Adaptable.
Screen Size	21.5"	21.5"	21.5"	15.6"
Probe Ports	4	3	3	1 (expandable to 3)
Touch Panel	●	●		
Height-Adjustable Control Panel	●	●		Optional with Advanced Cart
Swivel Control Panel	●			
Flexible monitor arm	●			
Gel Warmer	●	Optional		Optional with Advanced Cart
SSD	●	●	●	●
Windows® 10	●	●	●	●
Probes	18	12	9	12
Battery	Optional (15-Minute Capacity)	Optional (15-Minute Capacity)	Optional (15-Minute Capacity)	45 Minutes (More options available)
Wireless	Optional	●	●	Optional
DICOM®	●	●	●	●
Whizz Label	●	●		
Scan Coach	●	●	●	●
Scan Assistant	●	●	●	●
Whizz Dynamic Image Tuning	●	●	●	●
Whizz Easy Style	●			
Whizz Color Flow (Whizz CF)	●	●	●	
Follow-Up Tool	●	●	Optional	●
Needle Recognition	●	●	Optional	Optional
Auto Bladder	●	●	Optional	●
5-Year Warranty	●	●	●	●
4 Hours Remote eClinical Training	●	●	●	●
On-Site 1 Day Training	Optional	Optional	Optional	Optional
Flexible Financing Options	●	●	●	●

*Versana Premier V2, Versana Balance V2, Versana Essential R2 and Versana Active V1.5

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

©2023 GE HealthCare.

Versana Essential, Versana Premier, Versana Balance and Versana Active are trademarks of GE HealthCare. Windows is a registered trademark of Microsoft Corporation. DICOM is a registered trademark of the National Electrical Manufacturers Association. GE is a trademark of General Electric Company used under trademark license.

April 2023 JB01574XU

