

CRITIKON CLASSIC-CUF

Blood pressure cuffs

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong to withstand repeated inflations – yet economical enough for single patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE HealthCare patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling more accurate measurement

Constructed for patient safety and durability

- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment

Designed to improve workflow

- DINACLICK's easy-to-use design allows clinicians to quickly connect/disconnect adult and pediatric blood pressure cuffs
- Compliant with AAMI standard and AHA sizing guidelines
- DINACLICK is compatible with GE HealthCare and non-GE HealthCare automated BP devices, as well as on a variety of manual devices, including wall manometers and roll around sphygmomanometers, enabling a universal connection system throughout the facility

Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



Scan the QR code below to view our DINACLICK™ Connection System video



CLASSIC-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant
					
Size	Ink color	Limb circumference	2-Tube DINACLICK connector	1-Tube bayonet connector	2-Tube DINACLICK connector
Infant	Orange	8-13 cm	CLA-P1-2A-5S	2300	CLA-P1-2A
Child	Green	12-19 cm	CLA-P2-2A-5S	2301	CLA-P2-2A
Small adult	Light blue	17-25 cm	CLA-A1-2A-5S	2302	CLA-A1-2A
Adult	Navy	23-33 cm	CLA-A2-2A-5S	2303	CLA-A2-2A
Adult long	Navy	23-33 cm	CLA-A2-2A-L-5S	2306	CLA-A2-2A-L
Large adult	Rose	31-40 cm	CLA-A3-2A-5S	2304	CLA-A3-2A
Large adult long	Rose	31-40 cm	CLA-A3-2A-L-5S	2307	CLA-A3-2A-L
Thigh	Brown	38-50 cm	CLA-T1-2A-5S	2305	CLA-T1-2A

Note: 20 cuffs/box

DINACLICK assortment packs

Connector type	Part number (ISO 80369-5 compliant)	Part number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK connectors	2059304-501	2059304-001	5 small adult, 10 adult, 5 large adult (20 cuffs/box)
	2059304-502	2059304-002	1 small adult, 1 adult, 1 large adult (3 cuffs/pack)

To learn more about GE HealthCare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details. CLASSIC-CUF, CRITIKON, DINACLICK and Datex-Ohmeda are trademarks of GE HealthCare. GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

©2023 GE HealthCare.

January 2023
JB00412XX