

CRITIKON DURA-CUF

Blood pressure cuffs

CRITIKON™ DURA-CUF™ is gentle and designed for durability. Made with tough, flexible nylon, they are convenient and easy-to-clean with curved corners to enhance your patient's comfort. They are constructed for multiple patient use in low-acuity care settings and ideal for physicians offices where the risk of contamination is diminished.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large, bold print helps you quickly choose the right cuff, enabling a more accurate measurement

Constructed for patient safety and durability





- Soft, breathable, easy-to-clean material is not made with PVC, making them safer for the environment, or natural rubber latex, making them safer for your patients
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- DINACLICK easy-to-use design allows clinicians to quickly connect/disconnect adult and pediatric blood pressure cuffs
- Compliant with AAMI standard and AHA sizing guidelines
- DINACLICK is compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and roll around sphygmomanometers, enabling a universal connection system throughout the facility

Scan the QR code below to view our DINACLICK™ Connection System video



DURA-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant	
						
Size	Ink color	Limb circumference	2-Tube DINACLICK connector	1-Tube Bayonet connector	2-Tube DINACLICK connector	DURA inflation system cuffs
Infant	Orange	8-13 cm	DUR-P1-2A-5S	002274	DUR-P1-2A	002385
Child	Green	12-19 cm	DUR-P2-2A-5S	002275	DUR-P2-2A	002384
Child long	Green	12-19 cm	DUR-P2-2A-L-5S	—	DUR-P2-2A-L	002386
Small adult	Royal blue	17-25 cm	DUR-A1-2A-5S	002276	DUR-A1-2A	002383
Small adult long	Royal blue	17-25 cm	DUR A1-2A-L-5S	—	DUR A1-2A-L	002387
Adult w/hanger	Navy	23-33 cm	—	—	—	002380
Adult	Navy	23-33 cm	DUR-A2-2A-5S	002277	DUR-A2-2A	—
Adult long	Navy	23-33 cm	DUR-A2-2A-L-5S	002130	DUR-A2-2A-L	2388
Large adult	Wine	31-40 cm	DUR-A3-2A-5S	002278	DUR-A3-2A	002381
Large adult long	Wine	31-40 cm	DUR-A3-2A-L-5S	002131	DUR-A3-2A-L	2389
Thigh	Brown	38-50 cm	DUR-T1-2A-5S	002279	DUR-T1-2A	002382

Note: 5 cuffs/box

DINACLICK assortment packs

Connector Type	Part number (ISO 80369-5 compliant)	Part number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK connectors	2059301-501	2059301-001	1 small adult, 1 adult, 1 large adult (3 cuffs/box)
	2059301-502	2059301-002	2 infant, 3 child, 1 small adult (6 cuffs/box)
	2059301-503	2059301-003	1 small adult, 3 adult, 1 adult long, 1 large adult (6 cuffs/box)

To learn more about GE HealthCare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details. CRITIKON, DURA-CUF, DINACLICK and Datex-Ohmeda are trademarks of GE HealthCare. GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

©2023 GE HealthCare.

January 2023
JB00414XX

