

# CRITIKON SOFT-CUF

## Blood pressure cuffs

Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.



### Designed for patient safety, comfort and accuracy

- Priced economically for pairing with a single patient throughout their hospital stay to help with infection control efforts
- Soft, absorbent material is not made with DEHP or natural rubber latex to help enhance patient safety
- Designed with rounded corners for patient comfort
- Designed and clinically validated for performance accuracy to help you make well-informed decisions
- DINACLICK™ connector meets ISO 80369-5 international standards, lowering the risk of improper connections to other types of medical equipment
- A full range of cuff sizes are available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit

### Designed to improve workflow



- DINACLICK's easy-to-use design allows clinicians to quickly connect/disconnect adult and pediatric blood pressure cuffs
- Compliant with AAMI standard and AHA sizing guidelines
- DINACLICK is compatible with GE and non-GE automated BP devices, as well as on a variety of manual devices, including wall manometers and roll around sphygmomanometers, enabling a universal connection system throughout the facility

### Cuff assortment packs available

- Assortment packs are available to be flexible for all users and simplified ordering

Scan the QR code below to view our DINACLICK™ Connection System video



SOFT-CUF			ISO 80369-5 compliant		ISO 80369-1 compliant	
						
Size	Ink color	Limb circumference	2-Tube DINACLICK connector	1-Tube bayonet connector	2-Tube DINACLICK connector	SOFT inflation system cuffs
Infant	Orange	8-13 cm	SFT-P1-2A-5S	2360	SFT-P1-2A	—
Child	Green	12-19 cm	SFT-P2-2A-5S	2361	SFT-P2-2A	002584
Child long	Green	12-19 cm	SFT-P2-2A-L-5S	—	SFT-P2-2A-L	—
Small adult	Light blue	17-25 cm	SFT-A1-2A-5S	2362	SFT-A1-2A	002583
Small adult long	Light blue	17-25 cm	SFT-A1-2A-L-5S	—	SFT-A1-2A-L	—
Adult	Navy	23-33 cm	SFT-A2-2A-5S	2363	SFT-A2-2A	002580
Adult long	Navy	23-33 cm	SFT-A2-2A-L-5S	2140	SFT-A2-2A-L	002582
Large adult	Rose	31-40 cm	SFT-A3-2A-5S	2364	SFT-A3-2A	002581
Large adult long	Rose	31-40 cm	SFT-A3-2A-L-5S	2141	SFT-A3-2A-L	002585
Thigh	Brown	38-50 cm	SFT-T1-2A-5S	2365	SFT-T1-2A	—

Note: 20 cuffs/box

### DINACLICK assortment packs

Connector type	Part number (ISO 80369-5 compliant)	Part number (ISO 80369-1 compliant)	Description
2-Tube DINACLICK	2059303-501	2059303-001	5 small adult, 10 adult, 5 large adult (20 cuffs/box)
	2059303-503	2059303-003	3 small adult, 7 adult, 3 adult long, 3 large adult, 4 RADIAL-CUF (forearm) (20 cuffs/box)
	2059303-506	2059303-006	1 small adult, 1 adult, 1 large adult, 1 RADIAL-CUF (4 cuffs/pack)

To learn more about GE HealthCare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](https://cuffs.gehealthcare.com)

1. Infant BP Cuff for Datex-Ohmeda™ monitors has Male Slip Luer Connector.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details. CRITIKON, SOFT-CUF, DINACLICK and Datex-Ohmeda are trademarks of GE HealthCare. GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

©2023 GE HealthCare.

January 2023  
JB00418XX