






Look closer at a more secure, streamlined ultrasound *workflow*

The GE HealthCare Venue Family with Ready, Scan, Go reduces clinicians' administrative time* so they can focus on patient care, while administrators profit from seamless reliability, increased productivity, and accurate billing for insurance reimbursement.

More accuracy, more security, more efficiency, less time

-  **Easy logins**—Patient wristband scanning with barcode allows clinicians to remain focused on patient care instead of tedious and often distracting login process.
-  **No more ghost scans**—ViewPoint™ 6 integration helps Venue Family seamlessly track the clinician and patient journey for more reliable, more accurate billing, and optimal care.
-  **Better security and control**—Imprivata allows next-level flexibility to protect patient information.
-  **Reliable reimbursement**—Accuracy is everything when it comes to insurance reimbursements, so Ready, Scan, Go with the Venue Family keeps all of the key patient, clinician, and scanning data automatically integrated into each procedure.



Ready, Scan, Go™ is enabled with three key functions to enhance your workflow:

- Imprivata RFID Badge Reader
- Patient Barcode Scanner
- ViewPoint 6 POCUS Module—Streamlined Ultrasound Reporting

Imprivata, patient barcode scanner, and ViewPoint 6 can be used independently or together for a complete workflow solution. Using all three solutions together provides up to 79% time savings with connected, patient-centric collaboration.* This enables a more efficient workflow with the ability to store, report, and record POCUS studies directly from the ultrasound device.

The Venue Family




Venue™

Venue Go™

Venue Fit™

Accelerate ultrasound workflows at the point of care

A faster, more accurate, and more secure image is as easy as “Ready, Scan, Go”



Patients Ultrasound Exam

- Scan clinician ID badge
- Scan patient wristband
- Scan patient
- Populate POCUS report details on Venue
- Send Study to ViewPoint 6

ViewPoint 6

- On ViewPoint 6, finalize report
- ViewPoint 6 sends the report to EMR and image to PACS
- Perform QA

79% time-savings by implementing a total encounter-based workflow solution including RFID badge scanning for clinicians, barcode scanning of patient bracelet information and reporting*

Visit gehealthcare.com/readyscango to learn more about how you can start integrating it into your patients' care.



*Based on DOC2757861

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details. Ready, Scan, Go, Venue, Venue Go, Venue Fit and ViewPoint are trademarks of GE HealthCare. Imprivata is a registered trademark of Imprivata, Inc. in the U.S. and other countries. GE and the GE monogram are trademarks of General Electric Company. Used under trademark license. ©2023 GE HealthCare.